

ManaRay™

Red Light
and Near-Infrared
Therapy System



ManaRay Shoulder

 **MANAMED**

Recovery activated.

ManaRay™

LED System

Indications for Use

The ManaRay device is intended to emit energy in the near-infrared and visible red spectrum at wavelengths of 660nm, 850nm, and 940nm to provide topical heating for the purpose of elevating tissue temperature for the relief of muscle and joint pain, stiffness, arthritis pain, and muscle spasm; this increases local blood circulation, healing, and pain relief.



ManaRay Ankle



ManaRay Knee



ManaRay Universal



ManaRay Neck



ManaRay Hand



ManaRay Back

Product Description

The ManaRay device is a wearable LED-based system with adjustable hook and loop straps designed to deliver red and near-infrared light therapy. It includes:

- A control unit with USB-C charging port and power indicator lights.
- A flexible LED pad with 60 + 60 LEDs emitting at 660nm, 850nm, and 940nm

The device is powered by a built-in lithium-ion battery and is designed for short-duration therapy sessions at home or on-the-go. The ManaRay has a built-in timer for easy use.

Suggested HCPCS: E0221

Product Name	Part Number
ManaRay Ankle	MRAY-ANK001
ManaRay Back	MRAY-BK001
ManaRay Hand	MRAY-HND001
ManaRay Knee	MRAY-KN001
ManaRay Neck	MRAY-NK001
ManaRay Shoulder	MRAY-SH001
ManaRay Universal	MRAY-UNI001

For more information on suggested HCPCS and our PDAC letters, go to manamed.com/pdac-letters

2612 Sirius Drive, Denton, TX 76208

888-508-0712 (toll free)

support@Manamed.com www.Manamed.com

Manamed and the Manamed logo are trademarks of Manamed LLC.

TD-MRAY-SS Rev 1.0



Recovery activated