



Product Main Number 4118



Underpad, ABENA Abri-Soft Basic, 90x60cm, light blue, fluff/nonwoven/PE

- ✓ Absorbent fluff core
- ✓ Soft and smooth nonwoven surface
- ✓ Rapid liquid absorption
- ✓ Liquid safe back sheet and sealed edges for extra protection
- ✓ Available in different sizes



Product Description

ABENA Abri-Soft Basic is used as a protective, disposable underpad for e.g. beds and chairs. Suitable for incontinent users, but also for non-sterile procedures. The underpad has sealed edges for optimum leakage protection. The surface is covered with a soft nonwoven material for enhanced user comfort. The use of disposable underpads reduces the risk of contamination.

ABENA Abri-Soft Basic offers light to moderate absorbency and is available in different sizes.



Specifications

| | |
|----------------------------------------------------------------------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Base name | Underpad |
| Brand | ABENA |
| Sub-Brand | Abri-Soft Basic |
| Secondary color | Light blue |
| Features | Non-sterile |
| Single or multiple use | Single use |
| Material | PE, fluff, nonwoven |
| Length/depth | 90 cm |
| Width | 60 cm |
| Weight, net | 76 g |
| Tolerance | +/- 5 cm |
| Absorption Capacity | 1700 ml |
| Certificates | CE. FSC Mix. MD. Nordic Swan Ecolabel. UKCA. |
| CE Class (Medical Devices) | Class I |
| Product or test standards | ISO 11948-1 |
| Directives, regulations and acts | MDR (EU) 2017/745 |
| Because of the products composition it has an expiry date from the production date on: | 5 years |
| Storage Instructions | Store dry, room temperature and no direct sunlight. |
| Product Disposal Instructions | Can be disposed of with normal household waste sorted according to local regulations. If contamination is suspected, the product must be discarded as clinical waste. |
| Packaging Disposal Instructions | Can be recycled or incinerated. |

Instructions for use/application

The underpad must be placed on the bed or chair with the blue side down. Ensure that the underpad is fully smoothed out before use.

Packaging data

| Unit | Contains | Length | Width | Height | EAN |
|------|----------|---------|---------|---------|---------------|
| cII | 4 pck | 60.9 cm | 40.4 cm | 33.3 cm | 5703538933752 |
| pck | 30 pcs | 28 cm | 20 cm | 32 cm | 5703538933738 |



Nordic Swan Ecolabel is an official Scandinavian environmental label. Products with the Nordic Swan Ecolabel meet stringent environmental requirements with regard to raw material, production, use, documentation, and safety.



FSC™ (Forest Stewardship Council) is an international non-governmental organization dedicated to the promotion of responsible forest management worldwide. FSC™-certified products are made from raw materials from well-managed forests and responsible sources. The label's traceability and integrity are upheld through stringent control measures, and by choosing FSC™-certified products, you help take care of forests and the people and wildlife that call them home. FSC™ Mix means the product consists of a blend of FSC certified wood and recycled materials. For more information, go to www.fsc.org.



Regulation (EU) 2017/745 is a regulation of the European Union on the clinical investigation and sale of medical devices for human use. It repeals Directive 93/42/EEC, which concerns medical devices, and Directive 90/385/EEC, which concerns active implantable medical devices, on 26 May 2021.



The CE mark guarantees that a product is safe to use and complies with all safety precautions. CE stands for Conformité Européenne (European Conformity) and is mainly found on electronic equipment, safety equipment, construction products, and medical devices.



The UKCA mark guarantees that a product is safe to use and complies with all safety precautions in line with product standards and legislation of the UK.



Your guarantee that our products have been tested for ingredients which can potentially cause skin irritation. Tests conducted by an independent third-party laboratory.



The product is only for single-use.



The product has a very soft surface for a more comfortable fit.



Raw material that absorbs its own weight 10 times and has excellent fluid dispersion capabilities.

