



Item number 226805

Irrigation solution, 30 ml, sodium chloride, isotonic, sterile



Product Description

Sterile irrigation solution made from isotonic NaCl 0.9%. Used for irrigating wounds and eyes. For cleaning wounds and skin. Can also be used for irrigation of eyes, filling water evaporators, and lavage and cleaning of instruments. Do not use for IV use. Single use only. Discard after use. Do not heat in microwave oven.



Specifications

Base name	Irrigation solution
Features	Sodium chloride, isotonic, sterile
Single or multiple use	Single use
Ingredients / Composition	Sodium chloride 0.9%, PE bottle, sodium chloride 9g water for injection to 1000 ml. Osmolality approx./n.300 mosm/kg H ₂ O. pH approx./n.6. ISOTONIC.
Volume, net	30 ml
Certifications	CE
CE Class (Medical Devices)	Class IIa
Directives, regulations and acts	93/42/EEC
Safety Instructions and Warnings	Do not use for injection. Use only if the liquid is clear and the container intact.
Because of the products composition it has an expiry date from the production date on:	3 years
Shelf Life, opened	24 hours handled clean, for the same patient
Storage Instructions	Store dry, clean and at room temperature.
Product Disposal Instructions	Dispose with household waste.
Packaging Disposal Instructions	Can be recycled or incinerated.

Packaging data

Unit	Contains	Length	Width	Height	EAN
cil	240 pcs	32 cm	26 cm	28 cm	7053249369080
pck	20 pcs	15 cm	12.5 cm	9 cm	7053249369080



The CE mark guarantees that a product is safe to use and complies with all safety precautions. CE stands for Conformité Européenne (European Conformity) and is mainly found on electronic equipment, safety equipment, construction products, and medical devices.

